
Setting up Micro-manager - additional information

- 1. After building the controller connect it via USB to your computer that has Micro manager installed
- 2. Flash the Teensy with the provided script (Check PJRC.com for further instructions)
 - a. Version 5 is with physical resistors depicted in Figure 1B
 - b. Version 6 utilized the internal Arduino-Teensy pull-up resistors via code
- 3. Install Micro-manager
- 4. Start-up Micro-Manager
- 5. When prompted at start-up to select a hardware configuration choose "Create new configuration". Alternatively, if Micro-manager is already running open the *Hardware configuration wizard* and "Create new configuration".
- 6. Click Next and select Arduino hub and click done.





- 7. Add a camera that is attached to the system and name it accordingly.
- 8. Click on done and then select the "Arduino" and "scan ports". (Step 3/10 of configuration wizard) a number of values will be returned (Step 4/10).
- 9. Click Next. If properly connected properties of the Arduino will show up.
- 10. Activate the **Arduino Switch state** and **Arduino Shutter** by ticking the box.



11. Step 6-8 click next.


12. Step 9 of 10 you can define names for the Arduino Switch states.

Switch states:

Switch	Function	Proposed name
State		
0	μManager off	External
1	Control LED 1 with μManager – Continuous light	Bright
2	Control LED 2 with µManager – Continuous light	Fluorescence
3	Pulsed LED 1 control LED 2 with foot switch	Bright - Foot Fluo
4	Alternate LED1 and LED 2 at 10 Hz, activate LED 2 with foot switch ¹	Dual overlay (10 Hz ['])
5	Alternate LED1 and LED 2 at 20 Hz, activate LED 2 with foot switch ¹	Dual overlay (20 Hz ['])

Notes:

- *I.* State 4 and 5 require an additional bean shell script to be executed for properly assigning look up tables to the two channels.
- *II.* Other switch states are not assigned.
 - 13. Step 10 of 10: Save the configuration file and click finish.
 - **14.** In the Micromanager window add a new group by clicking on the plus button here called **Imaging.**




15. In the **Imaging** Group all settings are pre-defined that will be displayed in the Preset window. For other cameras than the Photometrics CoolSnap, Exposure, trigger and binning mode might be sufficient.

Activate Camera Properties (here CoolSnap):

- a. Binning
- b. Clear Mode (needed for the PVcam/Photometrics cameras).
- c. Exposure
- d. Trigger Mode

Activate Arduino Properties:

- e. Switch State and
- f. Arduino Shutter



16. Define the Pre-sets for a quick selection of Switch States – here exemplified for Switch state 2 – Fluorescence or Fluo (see image).



Running the Bean Shell script

- Switch State 4 and 5 require the execution of a bean shell script.
- Load the script for 10 or 20 Hz in the bean shell environment of Micro-manager and execute before starting "Live" in the Micro-manager control window.

References

Bosse, J. B., Tanneti, N. S., Hogue, I. B., and Enquist, L. W. (2015). Open LED Illuminator: A Simple and Inexpensive LED Illuminator for Fast Multicolor Particle Tracking in Neurons. PLoS ONE 10, e0143547. doi:10.1371/journal.pone.0143547.s008.